



Product Datasheet

RPTOR antibody (orb306096)

Description	Mouse monoclonal antibody to RPTOR
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	RPTOR
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	149038
Uniprot ID	Q8N122
Dilution Range	IHC-P: 1:25, WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa, MCF-7, mo...