



## Product Datasheet

### VAMP8 antibody (orb306098)

|                            |                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to VAMP8                                                                                                                                  |
| <b>Species/Host</b>        | Mouse                                                                                                                                                               |
| <b>Reactivity</b>          | Human                                                                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                        |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                           |
| <b>Target</b>              | VAMP8                                                                                                                                                               |
| <b>Form/Appearance</b>     | Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS. |
| <b>Note</b>                | For research use only                                                                                                                                               |
| <b>Isotype</b>             | IgG1,k                                                                                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                          |
| <b>MW</b>                  | 11438                                                                                                                                                               |
| <b>Uniprot ID</b>          | <b>Q9BV40</b>                                                                                                                                                       |
| <b>Dilution Range</b>      | IHC-P: 1:25, WB: 1:1000-1:2000                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                     |



Western blot analysis of HeLa, THP-1, A4...