

Product Datasheet

PRKAA1 antibody (orb306102)

Description	Rabbit polyclonal antibody to PRKAA1
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	PRKAA1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64009
Uniprot ID	Q13131
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

