



Product Datasheet

CHRNA3 antibody (orb306132)

Description	Rabbit polyclonal antibody to CHRNA3
Species/Host	Rabbit
Reactivity	Hamster
Conjugation	Unconjugated
Tested Applications	WB
Target	CHRNA3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52729
Uniprot ID	Q05901
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart, brain...