




Product Datasheet

CLN3 antibody (orb306142)

Description	Rabbit polyclonal antibody to CLN3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	CLN3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47623
Uniprot ID	Q13286
Dilution Range	IHC-P: 1:25, FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of A431, KG-1, SH-...