

Product Datasheet

FLT3 antibody (orb306160)

Description	Rabbit polyclonal antibody to FLT3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	FLT3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	112903
Uniprot ID	P36888
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

