

Product Datasheet

PPP2R2C antibody (orb306182)

Description	Rabbit polyclonal antibody to PPP2R2C
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	PPP2R2C
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51515
Uniprot ID	Q9Y2T4
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

