



Product Datasheet

PCGF6 antibody (orb306185)

Description	Rabbit polyclonal antibody to PCGF6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	PCGF6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39047
Uniprot ID	Q9BYE7
Dilution Range	IF/ICC: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of human kidney ti...