

Product Datasheet

COX11 antibody (orb306187)

Description	Rabbit polyclonal antibody to COX11
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	COX11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31430
Uniprot ID	Q9Y6N1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

