



Product Datasheet

PPP2R5B antibody (orb306194)

Description	Rabbit polyclonal antibody to PPP2R5B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	PPP2R5B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57393
Uniprot ID	Q15173
Dilution Range	IF/ICC: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of SK-N-MC cell li...