




Product Datasheet

CHRM2 antibody (orb306256)

Description	Mouse monoclonal antibody to CHRM2
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	CHRM2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	51715
Uniprot ID	P08172
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, WB: 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of paraffin-embedded tissue.