

Product Datasheet

HTR4 antibody (orb30626)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target HTR4

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 43761

Uniprot ID [Q13639](#)

NCBI [NP_001035263.1](#), [NP_955525.1](#), [NP_001035262.2](#), [NP_000861.1](#), [NP_001035259.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

