



Product Datasheet

CHMP3 antibody (orb306311)

Description	Rabbit polyclonal antibody to CHMP3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	CHMP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25073
Uniprot ID	Q9Y3E7
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain and...