

Product Datasheet

CHST11 antibody (orb306326)

Description	Rabbit polyclonal antibody to CHST11
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	CHST11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41555
Uniprot ID	Q9NPF2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

