

Product Datasheet

ATP5D antibody (orb306367)

Description	Rabbit polyclonal antibody to ATP5D
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	ATP5F1D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17490
Uniprot ID	P30049
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.

