

Product Datasheet

B4GALT7 antibody (orb306389)

Description	Rabbit polyclonal antibody to B4GALT7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	B4GALT7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37406
Uniprot ID	Q9UBV7
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

