




Product Datasheet

SNCA antibody (orb306394)

Description	Rabbit polyclonal antibody to SNCA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	SNCA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14460
Uniprot ID	P37840
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sk...



Flow cytometric analysis of SH-SY5Y cell...