



Product Datasheet

Noto antibody (orb306398)

Description	Rabbit polyclonal antibody to Noto
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	Noto
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26322
Uniprot ID	Q5TIS6
Dilution Range	WB: 1:4000, WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...