




Product Datasheet

FGF9 antibody (orb306400)

Description	Rabbit polyclonal antibody to FGF9
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Target	FGF9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23441
Uniprot ID	P31371
Dilution Range	FACS: 1:25, IHC-P1:25, WB: 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of All lanes using ...



Immunohistochemical staining of human li...