

Product Datasheet

ACAP2 antibody (orb306408)

Description	Rabbit polyclonal antibody to ACAP2
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	ACAP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88029
Uniprot ID	Q15057
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.

