

Product Datasheet

OR5M9 antibody (orb306414)

Description	Rabbit polyclonal antibody to OR5M9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	OR5M9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35093
Uniprot ID	Q8NGP3
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.

