



Product Datasheet

FLT1 antibody (orb306443)

Description	Mouse monoclonal antibody to FLT1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	FLT1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	150769
Uniprot ID	P17948
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of All lanes using ...