



Product Datasheet

Ptk7 antibody (orb306465)

Description	Rabbit polyclonal antibody to Ptk7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	Ptk7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	117532
Uniprot ID	Q8BKG3
Dilution Range	WB: 1:2000, WB: 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3 whole c...