



Product Datasheet

NFKBIE antibody (orb306489)

Description	Rabbit polyclonal antibody to NFKBIE
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	NFKBIE
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52864
Uniprot ID	000221
Dilution Range	WB: 1:2000, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of All lanes using ...