

Product Datasheet

RFX3 antibody (orb30649)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target RFX3

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 83530

Uniprot ID [P48380](#)

NCBI [NP_002910.1](#), [NP_001269045.1](#), [NP_001269046.1](#), [NP_602304.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

