




Product Datasheet

Sox17 antibody (orb306495)

Description	Rabbit polyclonal antibody to Sox17
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	Sox17
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44646
Uniprot ID	Q61473
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of All lanes using ...



Immunohistochemical staining of mouse te...