




## Product Datasheet

### Shh antibody (orb306498)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Shh                                                                                                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | FC, IF, IHC-P, WB                                                                                                                                                            |
| <b>Target</b>              | Shh                                                                                                                                                                          |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>MW</b>                  | 47773                                                                                                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q62226</a>                                                                                                                                                       |
| <b>Dilution Range</b>      | WB: 1:2000, IF/ICC: 1:25, IHC-P: 1:25, FACS: 1:25                                                                                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |



Western blot analysis  
of All lanes using ...



Immunohistochemical  
staining of mouse  
ce...