

Product Datasheet

fabF antibody (HRP) (orb30655)

Description	fabF antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Staphylococcus aureus 3-oxoacyl-[acyl-carrier-protein] synthase 2 protein (1-414AA)
Target	fabF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5HHA1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.