

Product Datasheet

CLEC11A antibody (orb306563)

Description	Rabbit polyclonal antibody to CLEC11A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	CLEC11A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35695
Uniprot ID	Q9Y240
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

