




Product Datasheet

Pou5f1 antibody (orb306605)

Description	Rabbit polyclonal antibody to Pou5f1
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	Pou5f1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38216
Uniprot ID	P20263
Dilution Range	IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of F9 whole cell l...



Immunohistochemical staining of mouse br...