

Product Datasheet

DNM1L antibody (orb306609)

Description	Rabbit polyclonal antibody to DNM1L
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	DNM1L
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	81877
Uniprot ID	000429
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

