



Product Datasheet

ITGA3 antibody (orb306623)

Description	Rabbit polyclonal antibody to ITGA3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	ITGA3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	116612
Uniprot ID	P26006
Dilution Range	WB: 1:2000, WB: 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of All lanes using ...