



Product Datasheet

NFE4 antibody (orb306646)

Description	Rabbit polyclonal antibody to NFE4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	NFE4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19019
Uniprot ID	Q86UQ8
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of human kidney ly...