



Product Datasheet

DYNC1H1 antibody (orb306656)

Description	Rabbit polyclonal antibody to DYNC1H1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	DYNC1H1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	532408
Uniprot ID	Q14204
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...