



Product Datasheet

FN3KRP antibody (orb306690)

Description	Rabbit polyclonal antibody to FN3KRP
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	FN3KRP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34412
Uniprot ID	Q9HA64
Dilution Range	IHC-P: 1:25, WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human St...