



Product Datasheet

RAPGEF1 antibody (orb306693)

Description	Rabbit polyclonal antibody to RAPGEF1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	RAPGEF1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	120548
Uniprot ID	Q13905
Dilution Range	WB: 1:2000, WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of All lanes using ...