

Product Datasheet

RAB11FIP5 antibody (orb306724)

Description	Rabbit polyclonal antibody to RAB11FIP5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	RAB11FIP5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70415
Uniprot ID	Q9BXF6
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

