

Product Datasheet

MED22 antibody (orb30691)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MED22
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22221
Uniprot ID	Q15528
NCBI	NP_852468.1 , NP_598395.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

