

Product Datasheet

RPP38 antibody (orb30765)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RPP38
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31834
Uniprot ID	P78345
NCBI	NP_006405.2 , NP_892117.1 , NP_001252530.1 , NP_001091059.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

