

## Product Datasheet

### PAGE2B antibody (FITC) (orb30767)

|                        |                                                                          |
|------------------------|--------------------------------------------------------------------------|
| <b>Description</b>     | PAGE2B antibody (FITC)                                                   |
| <b>Species/Host</b>    | Rabbit                                                                   |
| <b>Reactivity</b>      | Human                                                                    |
| <b>Conjugation</b>     | FITC                                                                     |
| <b>Immunogen</b>       | Recombinant Human Putative G antigen family E member 3 protein (1-111AA) |
| <b>Target</b>          | PAGE2B                                                                   |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol           |
| <b>Form/Appearance</b> | Liquid                                                                   |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.            |
| <b>Note</b>            | For research use only                                                    |
| <b>Isotype</b>         | IgG                                                                      |
| <b>Clonality</b>       | Polyclonal                                                               |
| <b>Uniprot ID</b>      | <a href="#">Q5JRK9</a>                                                   |
| <b>Dilution Range</b>  | ELISA                                                                    |
| <b>Expiration Date</b> | 12 months from date of receipt.                                          |