



Product Datasheet

CTAGE5 antibody (orb30772)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	This CTAGE5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 775-804 amino acids from the C-terminal region of human CTAGE5.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	O15320
NCBI	NP_001234917.1 , NP_976229.1 , NP_001234918.1 , NP_976231.1 , NP_005921.2 , NP_001234919.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB435 cell ...