

Product Datasheet

AFAP1 antibody (orb30773)

Description	Rabbit polyclonal antibody to AFAP1 protein encoded by this gene is a Src binding partner. May play ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AFAP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	80725
Uniprot ID	Q8N556
NCBI	NP_940997.1 , NP_001128119.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

