



Product Datasheet

ASPH antibody (orb30803)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ASPH
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	85863
Uniprot ID	Q12797
NCBI	NP_001158227.1 , NP_115856.1 , NP_001158225.1 , NP_001158223.1 , NP_115857.1 , NP_001158222.1 , NP_115855.1 , NP_001158228.1 , NP_004309.2 , NP_064549.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cell line
...