




## Product Datasheet

### SUV39H2 antibody (orb30845)

|                            |                                                                                                                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GST Tag belongs to the GST super-family, and Pi family. It is a widel...                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                    |
| <b>Reactivity</b>          | Human, Rat                                                                                                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                              |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                                 |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                                         |
| <b>Target</b>              | SUV39H2                                                                                                                                                                                   |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS. |
| <b>Note</b>                | For research use only                                                                                                                                                                     |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                |
| <b>MW</b>                  | 46682                                                                                                                                                                                     |
| <b>Uniprot ID</b>          | <a href="#">Q9H5I1</a>                                                                                                                                                                    |
| <b>NCBI</b>                | <a href="#">NP_001180354.1</a> , <a href="#">NP_078946.1</a> ,<br><a href="#">NP_001180355.1</a> , <a href="#">NP_001180356.1</a> ,<br><a href="#">NP_001180353.1</a>                     |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:50-100                                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                           |



Immunohistochemical analysis of formalin...



WB analysis of A431  
K562 CEM cell line  
a...