

Product Datasheet

SH3D19 antibody (Biotin) (orb30849)

Description	SH3D19 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human SH3 domain-containing protein 19 protein (468-767AA)
Target	SH3D19
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5HYK7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.