

Product Datasheet

c-Kit antibody (orb308734)

Description

c-Kit antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P

Applications
Immunogen

KLH conjugated synthetic peptide derived from human c-Kit. Please contact us for the exact immunogen sequence. The peptide is available as orb374888.

Target

c-Kit. The antibody detects endogenous levels of c-Kit

Form/Appearance

Liquid: 10 mM PBS, pH 7.4 with 10 mg/ml rAlbumin, 0.03% Proclin 300 and 25% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

Uniprot ID

P10721, **P05532**

NCBI

000213, **000222**

Entrez


3815

Dilution Range


IF/ICC: 1:100-500, IHC-P: 1:500

Expiration Date


12 months from date of receipt.



IF analysis of human 293t cells tissue u...



Immunofluorescence image of human u251 c...



Immunohistochemical staining of human OV...