

Product Datasheet

EGFR isoform a variant antibody (orb308736)

Description

EGFR isoform a variant antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human EGFR isoform a variant. Please contact us for the exact immunogen sequence. The peptide is available as orb374890.

Target

EGFR isoform a variant. The antibody detects endogenous levels of EGFR isoform a variant

Form/Appearance

Liquid: 10 mM PBS, pH 7.4 with 10 mg/ml rAlbumin, 0.03% Proclin 300 and 25% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW


120 kDa

Uniprot ID
[Q9EQ09](#), [P00533](#), [O70156](#)
NCBI
[005219](#), [005228](#), [31](#), [21](#)
Entrez
[1956](#)
Dilution Range


IHC-P: 1:100-500, WB: 1:50-300

Expiration Date


12 months from date of receipt.



Immunofluorescence image of pig large intestine tissue ...



IF image of pig large intestine tissue ...



IHC-P image of mouse skin tissue using a...