

Product Datasheet

ompC antibody (orb308739)

Description	ompC antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from E.coli ompC, Please contact us for the exact immunogen sequence. The peptide is available as orb374893.
Target	ompC. The antibody detects endogenous levels of ompC
Form/Appearance	Liquid: 10 mM PBS, pH 7.4 with 10 mg/ml rAlbumin, 0.03% Proclin 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	40 kDa
Uniprot ID	P06996
Dilution Range	IHC-P: 1:100-300, WB: 1:100-500
Expiration Date	12 months from date of receipt.