



Product Datasheet

CD206 antibody (orb308746)


Description	Mouse monoclonal antibody to CD206
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, ICC, IHC-Fr, IP, WB
Immunogen	Purified human mannose receptor
Target	CD206
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P22897
Entrez	4360
Dilution Range	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml. Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD206 positive dendr...



Flow cytometry surface staining pattern ...



Western blotting analysis of human CD206...