


Product Datasheet


CD134 antibody (PE) (orb308749)

Description Mouse monoclonal antibody against CD134 conjugated to


Reactivity	Human, Primate
Conjugation	PE
Tested Applications	FC
Immunogen	HTLV 1-transformed HUT-102 cells
Target	CD134
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P43489
Entrez	7293
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD134
positive
CD25 ...



Flow
cytometry
surface
staining
pattern ...